



Dressderm 5-Layers Silicone Foam Dressing

Dressderm 5-Layer Silicone Foam Dressing Offers Superior Absorption, Protection, And Comfort For Moderate To Heavy Exuding Wounds

Features:

- Combines Multiple Layers To Manage Moisture, Protect The Wound, And Promote Healing
- Uses Gentle Silicone Adhesive For Secure Attachment And Easy Removal Without Causing Trauma
- Flexible Design Conforms To The Shape Of The Wound, **Ensuring Comfort And A Secure Fit**
- Ideal For Prolonged Wear, Reducing The Need For **Frequent Dressing Changes**



 Suitable For Managing Moderate To Heavy Exuding Wounds, Helps Maintain A Moist Healing Environment, Reduces The Risk Of Infection, And Minimizes Skin Damage

DRESSDERM STERILE TO PORT OF THE STERILE TO PORT OF THE



Directions:

- Disinfect The Wound With Sterile Water
- Apply The Dressing With The Wound Contact Layer **Facing The Wound**
- Use A Secondary Dressing To Secure It In Place

Storage:

- Store Below 25°c
- Keep Away From Sunlight And Water

Shelf Life:

• 3 Years

Code	Size
DFS0505	5 cm x 5 cm
DFS1010	10 cm x 10 cm
DFS1020	10 cm x 20 cm
DFS1025	10 cm x 25 cm
DFS1030	10 cm x 30 cm
DFS1212SS	12.5 cm x 12.5 cm (Sacrum)
DFS1818SL	18 cm x 18 cm (Sacrum)

*Other Sizes Available On Request

